
Antibodies Product Datasheet

Product Name:	TF Antibody
Catalog Number:	B0873
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Coagulation Factor III. AA range:246-295
Concentration:	1 mg/ml
Gene Name:	F3
Protein Name:	Tissue factor
Other Name:	F3; Tissue factor; TF; Coagulation factor III; Thromboplastin; CD antigen CD142
Human Swiss Prot No.:	P13726
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TF Polyclonal Antibody detects endogenous levels of TF protein.
Storage and Stability:	-20°C/1 year