

Antibodies Product Datasheet

Product Name:	FoxO3a (Phospho-Ser318/321) Antibody
Catalog Number:	A9703
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human FoxO3a (Ser318 and 321)
Concentration:	1 mg/ml
Gene Name:	FOXO3 FKHL1 FOXO3A
Protein Name:	FoxO3a (Ser318/321)
Other Name:	Forkhead box protein O3 (AF6q21 protein) (Forkhead in rhabdomyosarcoma-like 1)
Human Swiss Prot No.:	O43524
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat FoxO3a (phospho-Ser318 or 321)
Storage and Stability:	-20°C/1 year