

Antibodies Product Datasheet

Product Name:	Neuro D2 Antibody
Catalog Number:	C11108
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NDF2. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	NEUROD2
Protein Name:	Neurogenic differentiation factor 2
Other Name:	NEUROD2; BHLHA1; NDRF; Neurogenic differentiation factor 2; NeuroD2; Class A basic helix-loop-helix protein 1; bHLHa1; NeuroD-related factor; NDRF
Human Swiss Prot No.:	Q15784
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Neuro D2 Polyclonal Antibody detects endogenous levels of Neuro D2 protein.
Storage and Stability:	-20°C/1 year