

Antibodies Product Datasheet

Product Name:	Anti-Human Cytokeratin 16 (ABT059) Mouse Monoclonal Antibody
Catalog Number:	V0295
Clone ID:	ABT059
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	KRT16
Gene Name:	P08779
Other Name:	CK 16, CK-16, CK16, Cytokeratin 16, Cytokeratin-16, Cytokeratin16, FNEPPK, Focal non epidermolytic palmoplantar keratoderma, K 16, K16, K1C16_HUMAN, K1CP, Keratin 1 type I, Keratin 16, Keratin, Keratin type I cytoskeletal 16, Keratin-16, Keratin16, KRT 16, Krt16, KRT16A, NEPPK, PC1, type I cytoskeletal 16
Isotype	IgG2b, Kappa
Localization:	Tonsil
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CK16 protein, and shows no cross reaction with CK1, 5, 6, 7, 8,10, 13, 14, 15, 19, 20.
Storage and Stability:	-20°C/1 year