
Antibodies Product Datasheet

Product Name:	CDC25B Antibody
Catalog Number:	R12-2597
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	CDC25B Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDC25B.
Concentration:	1 mg/ml
Gene Name:	CDC25B
Other Name:	CDC25HU2; CDC25M2; Dual specificity phosphatase Cdc25B; M-phase inducer phosphatase 2; MPIP2
Human Swiss Prot No.:	P30305
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	CDC25B Antibodyantibody detects endogenous levels of CDC25B.
Storage and Stability:	Stable at -20°C for at least 1 year.