

## Antibodies Product Datasheet

<b>Product Name:</b>	RBCK1 Antibody
<b>Catalog Number:</b>	C48499
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human RBCK1
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	RBCK1 C20orf18 RNF54 UBCE7IP3 XAP3 XAP4
<b>Protein Name:</b>	RanBP-type and C3HC4-type zinc finger-containing protein 1 (EC 6.3.2.-) (HBV-associated factor 4) (Heme-oxidized IRP2 ubiquitin ligase 1) (HOIL-1) (Hepatitis B virus X-associated protein 4) (RING finger protein 54) (Ubiquitin-conjugating enzyme 7-interact
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q9BYM8
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of RBCK1 at Human, Mouse,Rat
<b>Storage and Stability:</b>	Store at -20° C/1 year