

Antibodies Product Datasheet

Product Name:	CDC2 Monoclonal Antibody
Catalog Number:	N0140
Reactivity:	Human
Applications:	WB,IF/ICC,FCM,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of CDC2 expressed in E. Coli.
Concentration:	
Gene Name:	CDK1
Protein Name:	Cell division protein kinase 1
Other Name:	CDK1; CDC2; CDC28A; CDKN1; P34CDC2; Cyclin-dependent kinase 1; CDK1; Cell division control protein 2 homolog; Cell division protein kinase 1; p34 protein kinase
Human Swiss Prot No.:	P06493
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CDC2 Monoclonal Antibody detects endogenous levels of CDC2 protein.
Storage and Stability:	-20°C/1 year