
Antibodies Product Datasheet

Product Name:	KLHL3 Antibody
Catalog Number:	R12-2956
Reactivity:	Human
Applications:	WB, IHC, ELISA
Dilutions:	KLHL3 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KLHL3.
Concentration:	1 mg/ml
Gene Name:	KLHL3
Other Name:	Kelch-like protein 3; KIAA1129
Human Swiss Prot No.:	Q9UH77/Q8N4I8
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	KLHL3 Antibody antibody detects endogenous levels of KLHL3.
Storage and Stability:	Stable at -20°C for at least 1 year.