

Antibodies Product Datasheet

Product Name:	ATRIP (Phospho-Ser224) Antibody
Catalog Number:	A8295
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATRIP around the phosphorylation site of Ser224. AA range:190-239
Concentration:	1 mg/ml
Gene Name:	ATRIP
Protein Name:	ATR-interacting protein
Other Name:	ATRIP; AGS1; ATR-interacting protein; ATM and Rad3-related-interacting protein
Human Swiss Prot No.:	Q8WXE1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-ATRIP (S224) Polyclonal Antibody detects endogenous levels of ATRIP protein only when phosphorylated at S224.
Storage and Stability:	-20°C/1 year