

Antibodies Product Datasheet

Product Name:	PACS2 Antibody
Catalog Number:	C42335
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human PACS2
Concentration:	1 mg/ml
Gene Name:	PACS2 KIAA0602 PACS1L
Protein Name:	Phosphofurin acidic cluster sorting protein 2 (PACS-2) (PACS1-like protein)
Other Name:	
Human Swiss Prot No.:	Q86VP3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human PACS2
Storage and Stability:	-20°C/1 year