

Antibodies Product Datasheet

Product Name:	CTRP8 Antibody
Catalog Number:	C15205
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human CTRP8.
Concentration:	1 mg/ml
Gene Name:	C1QTNF8
Protein Name:	Complement C1q tumor necrosis factor-related protein 8
Other Name:	C1QTNF8; Complement C1q tumor necrosis factor-related protein 8; C1q/TNF-related protein 8; CTRP8
Human Swiss Prot No.:	P60827
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CTRP8 Polyclonal Antibody detects endogenous levels of CTRP8 protein.
Storage and Stability:	-20°C/1 year