

Antibodies Product Datasheet

Product Name:	EP3 Antibody
Catalog Number:	G098
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PE2R3. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	PTGER3
Protein Name:	Prostaglandin E2 receptor EP3 subtype
Other Name:	PTGER3; Prostaglandin E2 receptor EP3 subtype; PGE receptor EP3 subtype; PGE2 receptor EP3 subtype; PGE2-R; Prostanoid EP3 receptor P43115
Human Swiss Prot No.:	P43115
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EP3 Polyclonal Antibody detects endogenous levels of EP3 protein.
Storage and Stability:	-20°C/1 year