

Antibodies Product Datasheet

Product Name:	Anti-Human Progesterone Receptor (ABT005) Mouse Monoclonal Antibody
Catalog Number:	V0503
Clone ID:	ABT005
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	PGR
Gene Name:	P06401
Other Name:	NR3C3, Nuclear receptor subfamily 3 group C member 3, PGR, PR, PRA, PRB, PRGR_HUMAN, Progesterone receptor, Progestin receptor form A, Progestin receptor form B
Isotype	IgG1, Kappa
Localization:	Cervix/ Breast carcinoma tissue which weak positive staining
Positive Tissue Control:	Nuclear
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human PR protein, including PR-A and PR-B. In western blotting of Hela cell lysate, there shows a 99 kDa band detected by the antibody.
Storage and Stability:	-20°C/1 year