

Antibodies Product Datasheet

Product Name:	IRS-1 Antibody
Catalog Number:	B7120
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IRS-1. AA range:274-323
Concentration:	1 mg/ml
Gene Name:	IRS1
Protein Name:	Insulin receptor substrate 1
Other Name:	IRS1; Insulin receptor substrate 1; IRS-1
Human Swiss Prot No.:	P35568
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IRS-1 Polyclonal Antibody detects endogenous levels of IRS-1 protein.
Storage and Stability:	-20°C/1 year