

Antibodies Product Datasheet

Product Name:	CK16 Monoclonal Antibody(6F6)
Catalog Number:	N1321
Reactivity:	Human, Mouse, Rat
Applications:	IHC-P, IF/ICC
Dilutions:	WB 500-2000 1:200 IF 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CK16
Concentration:	
Gene Name:	KRT16
Protein Name:	Keratin type I cytoskeletal 16
Other Name:	KRT16; KRT16A; Keratin, type I cytoskeletal 16; Cytokeratin-16; CK-16; Keratin-16; K16
Human Swiss Prot No.:	P08779
Formulation:	PBS, pH 7.4, containing 0.5% BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CK16 proteins.
Storage and Stability:	-20°C/1 year