

Antibodies Product Datasheet

Product Name:	C/EBP ε (Phospho-Thr74) Antibody
Catalog Number:	A0827
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human C/EBP-epsilon around the phosphorylation site of Thr74. AA range:40-89
Concentration:	1 mg/ml
Gene Name:	CEBPE
Protein Name:	CCAAT/enhancer-binding protein epsilon
Other Name:	CEBPE; CCAAT/enhancer-binding protein epsilon; C/EBP epsilon
Human Swiss Prot No.:	Q15744
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-C/EBP ϵ (T74) Polyclonal Antibody detects endogenous levels of C/EBP ϵ protein only when phosphorylated at T74.
Storage and Stability:	-20°C/1 year