

Antibodies Product Datasheet

Product Name:	LMO3 Antibody
Catalog Number:	C11382
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LMO3. AA range:96-145
Concentration:	1 mg/ml
Gene Name:	LMO3
Protein Name:	LIM domain only protein 3
Other Name:	LMO3; RBTN3; RBTN2; RHOM3; LIM domain only protein 3; LMO-3; Neuronal-specific transcription factor DAT1; Rhombotin-3
Human Swiss Prot No.:	Q8TAP4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LMO3 Polyclonal Antibody detects endogenous levels of LMO3 protein.
Storage and Stability:	-20°C/1 year