

Antibodies Product Datasheet

Product Name:	Caldesmon Antibody
Catalog Number:	C30874
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CALD1. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	CALD1
Protein Name:	Caldesmon
Other Name:	Caldesmon (CDM)
Human Swiss Prot No.:	Q05682
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Caldesmon Polyclonal Antibody detects endogenous levels of Caldesmon
Storage and Stability:	-20°C/1 year