

Antibodies Product Datasheet

Product Name:	Cleaved-ITI-H2 (D702) Antibody
Catalog Number:	L0300
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ITIH2. AA range:653-702
Concentration:	1 mg/ml
Gene Name:	ITIH2
Protein Name:	Inter-alpha-trypsin inhibitor heavy chain H2
Other Name:	ITIH2; IGHEP2; Inter-alpha-trypsin inhibitor heavy chain H2; ITI heavy chain H2; ITI-HC2; Inter-alpha-inhibitor heavy chain 2; Inter-alpha-trypsin inhibitor complex component II; Serum-derived hyaluronan-associated protein; SHAP
Human Swiss Prot No.:	P19823
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-ITI-H2 (D702) Polyclonal Antibody detects endogenous levels of fragment of activated ITI-H2 protein resulting from cleavage adjacent to D702.
Storage and Stability:	-20°C/1 year