
Antibodies Product Datasheet

Product Name:	Raf-1 Antibody
Catalog Number:	B7207
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human C-RAF. AA range:305-354
Concentration:	1 mg/ml
Gene Name:	RAF1
Protein Name:	RAF proto-oncogene serine/threonine-protein kinase
Other Name:	RAF1; RAF; RAF proto-oncogene serine/threonine-protein kinase; Proto-oncogene c-RAF; cRaf; Raf-1
Human Swiss Prot No.:	P04049
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Raf-1 Polyclonal Antibody detects endogenous levels of Raf-1 protein.
Storage and Stability:	-20°C/1 year