

Antibodies Product Datasheet

Product Name:	GRB2 mouse mAb
Catalog Number:	N0861
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant protein of human GRB2.
Concentration:	1 mg/ml
Gene Name:	grb2
Protein Name:	
Other Name:	Abundant SRC homology; Adapter protein GRB2; ASH; Ash protein; EGFRBP GRB2; Epidermal growth factor receptor binding protein; Epidermal growth factor receptor binding protein GRB2; GRB 2; GRB2 adapter protein; Grb2; GRB2_HUMAN; Grb3 3; Growth factor recep
Human Swiss Prot No.:	P62993
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GRB2 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year