

Antibodies Product Datasheet

Product Name:	RHOC Antibody
Catalog Number:	C48081
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RHOC
Concentration:	1 mg/ml
Gene Name:	RHOC ARH9 ARHC
Protein Name:	Rho-related GTP-binding protein RhoC (Rho cDNA clone 9) (h9)
Other Name:	
Human Swiss Prot No.:	P08134
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RHOC at Human, Mouse
Storage and Stability:	Store at -20° C/1 year