

Antibodies Product Datasheet

Product Name:	SLC40A1 Antibody
Catalog Number:	C51057
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SLC40A1. AA range 1-100
Concentration:	1 mg/ml
Gene Name:	SLC40A1 FPN1 IREG1 SLC11A3 MSTP079
Protein Name:	Solute carrier family 40 member 1 (Ferroportin-1) (Iron-regulated transporter 1)
Other Name:	
Human Swiss Prot No.:	Q9NP59
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SLC40A1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year