

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-p40 [ABT-P40]
<b>Catalog Number:</b>	V0144
<b>Clone ID:</b>	ABT-P40
<b>Reactivity:</b>	Human
<b>Applications:</b>	IHC
<b>Application Dilutions</b>	IHC 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	Q9H3D4
<b>Gene Name:</b>	TP63
<b>Other Name:</b>	KET, P63, P73H, P73L, TP73L
<b>Isotype</b>	Mouse IgG
<b>Localization:</b>	N/A
<b>Positive Tissue Control:</b>	N/A
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	Anti-p40 [ABT-P40] detects endogenous levels of total p40 protein.
<b>Storage and Stability:</b>	-20°C/1 year