
Antibodies Product Datasheet

Product Name:	CD86 Monoclonal Antibody
Catalog Number:	N0136
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human CD86 expressed in E. Coli.
Concentration:	
Gene Name:	CD86
Protein Name:	T-lymphocyte activation antigen CD86
Other Name:	CD86; CD28LG2; T-lymphocyte activation antigen CD86; Activation B7-2 antigen; B70; BU63; CTLA-4 counter-receptor B7.2; FUN-1; CD antigen CD86
Human Swiss Prot No.:	P42081
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	CD86 Monoclonal Antibody detects endogenous levels of CD86 protein.
Storage and Stability:	-20°C/1 year