
Antibodies Product Datasheet

Product Name:	CK17 Monoclonal Antibody(10A1)
Catalog Number:	N1346
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCIP
Dilutions:	WB: 1:1000 IP:1:200 IF 1:200 IHC 1:50-300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CK17
Concentration:	
Gene Name:	KRT17
Protein Name:	Keratin type I cytoskeletal 17
Other Name:	KRT17; Keratin, type I cytoskeletal 17; 39.1; Cytokeratin-17; CK-17; Keratin-17; K17
Human Swiss Prot No.:	Q04695
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects CK17 endogenous proteins.
Storage and Stability:	-20°C/1 year