

Antibodies Product Datasheet

Product Name:	DRS-1 Antibody
Catalog Number:	C17621
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PECL. AA range:73-122
Concentration:	1 mg/ml
Gene Name:	ECI2
Protein Name:	Enoyl-CoA delta isomerase 2 mitochondrial
Other Name:	ECI2; DRS1; HCA88; PECL; Enoyl-CoA delta isomerase 2; mitochondrial; DRS-1; Delta(3),delta(2)-enoyl-CoA isomerase; D3,D2-enoyl-CoA isomerase; Diazepam-binding inhibitor-related protein 1; DBI-related protein 1; Dodecenoyl-CoA isomerase; Hep
Human Swiss Prot No.:	O75521
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DRS-1 Polyclonal Antibody detects endogenous levels of DRS-1 protein.
Storage and Stability:	-20°C/1 year