

Antibodies Product Datasheet

Product Name:	B3GNT7 Antibody
Catalog Number:	C50710
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human B3GNT7
Concentration:	1 mg/ml
Gene Name:	B3GNT7
Protein Name:	UDP-GlcNAc:betaGal beta-1,3-N-acetylglucosaminyltransferase 7 (BGnT-7) (Beta-1,3-Gn-T7) (Beta-1,3-N-acetylglucosaminyltransferase 7) (Beta3Gn-T7) (EC 2.4.1.-)
Other Name:	
Human Swiss Prot No.:	Q8NFL0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of B3GNT7 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year