

Antibodies Product Datasheet

Product Name:	FKBP14 Antibody
Catalog Number:	C49089
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FKBP14
Concentration:	1 mg/ml
Gene Name:	FKBP14 FKBP22 UNQ322/PRO381
Protein Name:	Peptidyl-prolyl cis-trans isomerase FKBP14 (PPIase FKBP14) (EC 5.2.1.8) (22 kDa FK506-binding protein) (22 kDa FKBP) (FKBP-22) (FK506-binding protein 14) (FKBP-14) (Rotamase)
Other Name:	
Human Swiss Prot No.:	Q9NWM8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FKBP14 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year