

Antibodies Product Datasheet

Product Name:	MAST3 Antibody
Catalog Number:	C10132
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAST3. AA range:1031-1080
Concentration:	1 mg/ml
Gene Name:	MAST3
Protein Name:	Microtubule-associated serine/threonine-protein kinase 3
Other Name:	MAST3; KIAA0561; Microtubule-associated serine/threonine-protein kinase 3
Human Swiss Prot No.:	O60307
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MAST3 Polyclonal Antibody detects endogenous levels of MAST3 protein.
Storage and Stability:	-20°C/1 year