

Antibodies Product Datasheet

Product Name:	NFAT5 Mouse mAb(1D6)
Catalog Number:	N1453
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NFAT5
Concentration:	1 mg/ml
Gene Name:	NFAT5 KIAA0827 TONEBP
Protein Name:	NFAT5
Other Name:	Nuclear factor of activated T-cells 5 (NF-AT5) (T-cell transcription factor NFAT5) (Tonicity-responsive enhancer-binding protein) (TonE-binding protein) (TonEBP)
Human Swiss Prot No.:	O94916
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.26% sodium azide.
Specificity:	This antibody detects endogenous levels of NFAT5 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year