

Antibodies Product Datasheet

Product Name:	CLINT1 Antibody
Catalog Number:	C50881
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CLINT1
Concentration:	1 mg/ml
Gene Name:	CLINT1 ENTH EPN4 EPNR KIAA0171
Protein Name:	Clathrin interactor 1 (Clathrin-interacting protein localized in the trans-Golgi region) (Clint) (Enthoprotin) (Epsin-4) (Epsin-related protein) (EpsinR)
Other Name:	
Human Swiss Prot No.:	Q14677
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CLINT1 at Human
Storage and Stability:	Store at -20° C/1 year