
Antibodies Product Datasheet

Product Name:	α -tubulin Monoclonal Antibody
Catalog Number:	N1380
Reactivity:	Zebrafish
Applications:	WB
Dilutions:	WB: 1:500-10000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Tubulin alpha-1A chain
Concentration:	
Gene Name:	TUBA1A
Protein Name:	Tubulin alpha-1A chain
Other Name:	TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-III
Human Swiss Prot No.:	Q13509
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects Zebrafish endogenous α -tubulin protein.
Storage and Stability:	-20°C/1 year