
Antibodies Product Datasheet

Product Name:	Cytokeratin 15 Antibody
Catalog Number:	C0240
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Keratin 15. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	KRT15
Protein Name:	Keratin type I cytoskeletal 15
Other Name:	KRT15; KRTB; Keratin; type I cytoskeletal 15; Cytokeratin-15; CK-15; Keratin-15; K15
Human Swiss Prot No.:	P19012
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cytokeratin 15 Polyclonal Antibody detects endogenous levels of Cytokeratin 15 protein.
Storage and Stability:	-20°C/1 year