

Antibodies Product Datasheet

Product Name:	SLC37A1 Antibody
Catalog Number:	C50349
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SLC37A1
Concentration:	1 mg/ml
Gene Name:	SLC37A1 G3PP
Protein Name:	Glycerol-3-phosphate transporter (G-3-P transporter) (Glycerol-3-phosphate permease) (G-3-P permease) (Solute carrier family 37 member 1)
Other Name:	
Human Swiss Prot No.:	P57057
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SLC37A1 at Human
Storage and Stability:	Store at -20° C/1 year