

Antibodies Product Datasheet

Product Name:	PRMT4/CARM1 mouse mAb
Catalog Number:	N0780
Reactivity:	Human, Mouse
Applications:	WB, IP
Dilutions:	wb 1:200-1:500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human PRMT4/CARM1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	carm1
Protein Name:	
Other Name:	carm1; CARM1_HUMAN; Coactivator associated arginine methyltransferase 1; Coactivator-associated arginine methyltransferase 1; Histone arginine methyltransferase CARM 1; Histone arginine methyltransferase CARM1; Histone-arginine methyltransferase CARM1; PRMT 4; P
Human Swiss Prot No.:	Q86X55
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PRMT4/CARM1 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year