

Antibodies Product Datasheet

Product Name:	CtIP Antibody
Catalog Number:	B0416
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CTIP. AA range:293-342
Concentration:	1 mg/ml
Gene Name:	RBBP8
Protein Name:	DNA endonuclease RBBP8
Other Name:	RBBP8; CTIP; DNA endonuclease RBBP8; CtBP-interacting protein; CtIP; Retinoblastoma-binding protein 8; RBBP-8; Retinoblastoma-interacting protein and myosin-like; RIM; Sporulation in the absence of SPO11 protein 2 homolog; SAE2
Human Swiss Prot No.:	Q99708
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CtIP Polyclonal Antibody detects endogenous levels of CtIP protein.
Storage and Stability:	-20°C/1 year