

Antibodies Product Datasheet

Product Name:	PIP5K Antibody
Catalog Number:	R12-3327
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	PIP5K Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PIP5K.
Concentration:	1 mg/ml
Gene Name:	PIKFYVE
Other Name:	1- phosphatidylinositol-4-phosphate 5-kinase; FYV1; FYVE finger-containing phosphoinositide kinase; KIAA0981; PIKfyve; PIP5K3; PtdIns(4)P-5- kinase; p235
Human Swiss Prot No.:	Q9Y2I7
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	PIP5K Antibodyantibody detects endogenous levels of PIP5K.
Storage and Stability:	Stable at -20°C for at least 1 year.