

Antibodies Product Datasheet

Product Name:	MEMO1 Antibody
Catalog Number:	C50345
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MEMO1
Concentration:	1 mg/ml
Gene Name:	MEMO1 C2orf4 MEMO NS5ATP7 CGI-27
Protein Name:	Protein MEMO1 (C21orf19-like protein) (Hepatitis C virus NS5A-transactivated protein 7) (HCV NS5A-transactivated protein 7) (Mediator of ErbB2-driven cell motility 1) (Mediator of cell motility 1) (Memo-1)
Other Name:	
Human Swiss Prot No.:	Q9Y316
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MEMO1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year