

Antibodies Product Datasheet

Product Name:	DDX3 Antibody
Catalog Number:	R12-2111
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DDX3. AA range:14-63
Concentration:	1 mg/ml
Gene Name:	DDX3X/DDX3Y
Protein Name:	ATP-dependent RNA helicase DDX3X/Y
Other Name:	DDX3X; DBX; DDX3; ATP-dependent RNA helicase DDX3X; DEAD box protein 3; X-chromosomal; DEAD box, X isoform; Helicase-like protein 2; HLP2; DDX3Y; DBY; ATP-dependent RNA helicase DDX3Y; DEAD box protein 3, Y-chromosomal
Human Swiss Prot No.:	O00571/O15523
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DDX3 Polyclonal Antibody detects endogenous levels of DDX3 protein.
Storage and Stability:	-20°C/1 year