

Antibodies Product Datasheet

Product Name:	SENP8 Antibody
Catalog Number:	C0371
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SENP8. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	SENP8
Protein Name:	Sentrin-specific protease 8
Other Name:	SENP8; DEN1; NEDP1; PRSC2; FKSG8; Sentrin-specific protease 8; Deneddylase-1; NEDD8-specific protease 1; Protease; cysteine 2; Sentrin/SUMO-specific protease SENP8
Human Swiss Prot No.:	Q96LD8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SENP8 Polyclonal Antibody detects endogenous levels of SENP8 protein.
Storage and Stability:	-20°C/1 year