

Antibodies Product Datasheet

Product Name:	IGFBP3 (Phospho-Ser183) Antibody
Catalog Number:	A1052
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IGFBP-3 around the phosphorylation site of Ser183. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	IGFBP3
Protein Name:	Insulin-like growth factor-binding protein 3
Other Name:	IGFBP3; IBP3; Insulin-like growth factor-binding protein 3; IBP-3; IGF-binding protein 3; IGFBP-3
Human Swiss Prot No.:	P17936
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-IGFBP3 (S183) Polyclonal Antibody detects endogenous levels of IGFBP3 protein only when phosphorylated at S183.
Storage and Stability:	-20°C/1 year