

Antibodies Product Datasheet

Product Name:	CKR-5 Antibody
Catalog Number:	B0062
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CCR5. AA range:303-352
Concentration:	1 mg/ml
Gene Name:	CCR5
Protein Name:	C-C chemokine receptor type 5
Other Name:	CCR5; CMKBR5; C-C chemokine receptor type 5; C-C CKR-5; CC-CKR-5; CCR-5; CCR5; CHEMR13; HIV-1 fusion coreceptor; CD antigen CD195
Human Swiss Prot No.:	P51681
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CKR-5 Polyclonal Antibody detects endogenous levels of CKR-5 protein.
Storage and Stability:	-20°C/1 year