
Antibodies Product Datasheet

Product Name:	Rabbit Anti Rat IgG (H&L) - Alexa Fluor 680
Catalog Number:	SA1027
Reactivity:	Rat
Conjugate:	Alexa Fluor 680
Applications:	IF/ICC,FCM,WB
Dilutions:	Suggested starting 1:50-1:1000 dilutions for most fluorescent applications, while 1:5000-20000 dilutions for fluorescent assays with LI-1:COR® Odyssey equipment would be recommended.
Modification:	Not Modification Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Concentration:	1mg/mL
Other Name:	Rabbit Anti-Rat IgG (H&L) - Alexa Fluor 680
Formulation:	1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing