

Antibodies Product Datasheet

Product Name:	DDX56 Antibody
Catalog Number:	C50573
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DDX56
Concentration:	1 mg/ml
Gene Name:	DDX56 DDX21 NOH61
Protein Name:	Probable ATP-dependent RNA helicase DDX56 (EC 3.6.4.13) (ATP-dependent 61 kDa nucleolar RNA helicase) (DEAD box protein 21) (DEAD box protein 56)
Other Name:	
Human Swiss Prot No.:	Q9NY93
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DDX56 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year