

Antibodies Product Datasheet

Product Name:	ME3 Antibody
Catalog Number:	C16864
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ME3. AA range:545-594
Concentration:	1 mg/ml
Gene Name:	ME3
Protein Name:	NADP-dependent malic enzyme mitochondrial
Other Name:	ME3; NADP-dependent malic enzyme; mitochondrial; NADP-ME; Malic enzyme 3
Human Swiss Prot No.:	Q16798
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ME3 Polyclonal Antibody detects endogenous levels of ME3 protein.
Storage and Stability:	-20°C/1 year