

Antibodies Product Datasheet

Product Name:	TudorSN Antibody
Catalog Number:	C20602
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human TudorSN.
Concentration:	1 mg/ml
Gene Name:	SND1
Protein Name:	Staphylococcal nuclease domain-containing protein 1
Other Name:	SND1; TDRD11; Staphylococcal nuclease domain-containing protein 1; 100 kDa coactivator; EBNA2 coactivator p100; Tudor domain-containing protein 11; p100 co-activator
Human Swiss Prot No.:	Q7KZF4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TudorSN Polyclonal Antibody detects endogenous levels of TudorSN protein.
Storage and Stability:	-20°C/1 year