

Antibodies Product Datasheet

Product Name:	HLA-DPB1 Antibody
Catalog Number:	C50834
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HLA-DPB1
Concentration:	1 mg/ml
Gene Name:	HLA-DPB1 HLA-DP1B
Protein Name:	HLA class II histocompatibility antigen, DP beta 1 chain (HLA class II histocompatibility antigen, DP(W4) beta chain) (MHC class II antigen DPB1)
Other Name:	
Human Swiss Prot No.:	P04440
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HLA-DPB1 at Human
Storage and Stability:	Store at -20° C/1 year