

Antibodies Product Datasheet

Product Name:	IL-1F5 Antibody
Catalog Number:	C30617
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IL36RN. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	IL36RN
Protein Name:	Interleukin-36 receptor antagonist protein
Other Name:	IL36RN; FIL1D; IL1F5; IL1HY1; IL1L1; IL1RP3; Interleukin-36 receptor antagonist protein; FIL1 delta; IL-1-related protein 3; IL-1RP3; Interleukin-1 HY1; IL-1HY1; Interleukin-1 delta; IL-1 delta; Interleukin-1 family member 5; IL-1F5; Interleukin-1 recepto
Human Swiss Prot No.:	Q9UBH0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-1F5 Polyclonal Antibody detects endogenous levels of IL-1F5 protein.
Storage and Stability:	-20°C/1 year