
Antibodies Product Datasheet

Product Name:	P70 S6 Kinase mouse Monoclonal Antibody(1C7)
Catalog Number:	N1586
Reactivity:	Human,Rat,Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of P70 S6 Kinase of RPS6KB1
Concentration:	1 mg/ml
Gene Name:	RPS6KB1
Protein Name:	RPS6KB1
Other Name:	RPS6KB1
Human Swiss Prot No.:	P23443
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	P70 S6 Kinase protein detects endogenous levels of RPS6KB1
Storage and Stability:	-20°C/1 year